Stereoscopic Micro Particle Image Velocimetry

Technische Universität Braunschweig | Institute for Chemical und Thermal Process Engineering ictv@tu-braunschweig.de | Phone +49 (0) 531 391-2791

Principle of Measuring

the tracers are excited to glow.

appears as bright dots on the images.

of the tracers.

calculated.

Time between pulses

Trigger rate

■ The laser (Nd:YAG) fires a laser pulse for a few nanoseconds with a

defined wavelength of 532 nm, adapted to the excitation of fluorescence

A dichroic mirror reflects the laser beam into the flow channel in which

The light emitted by the tracers has a different wavelength (584 nm) than the laser and can transmit through the dichroic mirror.

The relevant parameters, i.e., the time between two laser pulses, the

With the divergent position of the particles on each image of the image

trigger rate and the number of images must be defined.

pair and the time between the images the velocity vectors are

The emitted light it is detected by the image sensors of the cameras and

Objectives

- Three-dimensional visualization of streamlines (2D3C).
- Quantification of flow vectors within the velocity fields.
- Determination of local turbulence intensities.
- Calculation of wall shear stress.





Image Processing

- Background subtraction to reduce wall reflection is applied (mean value filter).
- Analysis with adaptive correlation, which iteratively adjusts the shape and size of interrogation areas.
- Merged image are created from the calculated vectors of the single images.

Requirements

- Field of view:
- $50 \times 50 \ \mu\text{m}^2$ to $18 \times 18 \ \text{mm}^2$ Depth of field: up to 1 μm
- Depth of held, up to 1 pm
 Duration: ~ 60 ms to a few minutes
- Optical accessibility is required
- Tracers must be dispersed in the fluid (d_p = 1 20 μm, Ψ≈ 1)



Time between pulses

Number of images





Technische Universität Braunschweig

